


User Instruction for Gap Front Shower Stool S11



1. Introduction

This product is designed to be used to help aid showering; featuring a reinforced framework that caters for weight up to 154 kg. The clip-on perforated seat and the chair's integral handles enhance the users support whilst showering. The legs are height adjustable and the non-slip rubber feet ensure frame and patient stability.

Product Dimensions	
Seat height (from-to)	490mm-641mm
Seat width	440mm
Seat depth	370mm
Seat to top of arms	150mm
Width between arms	562mm
Footprint	380mm x 560mm

2. Contents of the Box

- Gap Front Shower Stool (S11)
- Gap Front Seat (72030H/L)
- Extended Legs & Pin Clips
- User Instructions (this leaflet)

3. Assembly Instructions


- Carefully unpack and examine the product to ensure it is complete.
- Before use check for any signs of damage e.g. cracked weld, flaking of paint, bent tube. If any signs occur please contact your supplier to exchange and refrain from using the product. **DO NOT ATTEMPT TO USE IF PRODUCT IS DAMAGED.**

4. Correct Use

- Drainage holes comply with MHRA guidelines.

Adjustable Height:

- Remove the pin clips from legs.



User Instruction for Gap Front Shower Stool S11

- Select the required seat height by gently sliding out each extension leg from the leg of the stool frame.
- At the desired height refit the pin clip. Ensure the pin clip fully goes through each leg and its extension leg.
- The other legs are adjusted in the same manner. It is important to ensure all legs are set to the same height.
- The stool is now ready for use.

5. Cleaning

- Use mild detergent with warm water to wipe down framework every time you use.
- Do not use an abrasive cleaner as this will damage the frame and nullify your warranty.

Use of power/jet washing & industrial washing machines:

- If either of the above is used to clean equipment care must be taken since high-pressure water can be penetrative.
- Care should be taken that frames are not placed too close to high-pressure water jets.
- After cleaning ensure any residual water is removed from the framework as it could cause permanent rusting of the frame (leave frames open in a warm environment for several hours to ensure all internal moisture evaporates).
- We cannot offer advice on industrial washing machines with regards to specific machine settings, as these will be dependant on the manufacturer's instructions.

6. Maintenance Requirements

The product should be checked on a regular basis for damage or wear and tear – this includes flaking or cracking paint. Any undue signs of damage should be noted and reported to the supplier in the first instance. Should damage be evident, then the user should refrain from using the product until advice has been sought from the supplier.

7. Warning

- This product is intended for users of a maximum of 154kg/25stone in weight.
- We recommend a full risk assessment should be carried out before use.

8. Warranty

This product carries a 12 month warranty from the date of purchase against manufacturing defects. Misuse or damage is not covered. This warranty does not affect your statutory rights. As manufacturers we cannot specify expected product life as this is based on the level of usage.

9. Disposal Instructions

- Appropriate hygiene standards should be applied before disposal to reduce the risk of infection.
- Materials used are essentially non hazardous and can be disposed of appropriately.

10. Contact Details

We have attempted to cover most eventualities, however if you have any specific queries, please contact:

Cefndy Healthcare

Cefndy Road

Rhyl

Denbighshire

LL18 2HG

Tel: 0044 1745 343877 (Quality Assurance Dept)

E-mail: cefndy.sales@denbighshire.gov.uk